Platform	Models	Picture	Material/Ctg	Asph/n/D/A	Notes			
Monofocal								
Tecnis (Tecnis, J&J)	ZCB00 DCB00 (preload) ZCUx/ZCTx (toric)	Q	Acrylic 1-Piece (Bag)  ©/©	-0.27 / 1.47 +5.0 to +34.0 119.3/118.8	Diameter is 6.0 mm but for powers > +21.00 D the optic is 4.9 mm ("U" for unpolished haptic) -0.27 µm of spherical aberration (SA)			
AcrySof (Alcon)	SA60WF/SA60AT SN60WF/SN60AT SA6AT <i>x</i> (toric) SN6AT <i>x</i> (toric)		Acrylic 1-Piece (Bag)  ©/©	-0.20/sph/1.55 +6.0 to +34.0 -AT to +40.0 119.0/118.7	Yellow (Natural) tint blocks blue light, Will be phased out by Alcon for Clareon "AT" series is spherical, <sup>(B)</sup> is used for +35 – +40 "WF" is aspheric			
Clareon (Alcon)	CC60WF/SY60WF CCW0Tx/CNW0Tx CAA0T0/CNA0T0		Acrylic 1-Piece (Bag) ©/©	-0.20 / 1.55 +6.0 to +30.0 119.1/118.8	Replacing the AcrySof platform; 6.0 mm optic -A0T0 as a CO <sub>2</sub> powered injector (AutonoMe) for one hand use; "glistening-free" material			
enVista (B+L)	MX60E MX60PL <sub>(preload)</sub> MX60ET <sub>(toric)</sub> MX60PT <sub>(preload toric)</sub>		Acrylic 1-Piece (Bag) ©/©	NEUTRAL/1.53 0.0 to +34.0 119.1/118.7	Doesn't change the asphericity of the eye's optical system; has eyelets for scleral fixation; large range of spherical diopteric powers Toric powers down to 1.25D; 6.0 mm optic			
Light Adjustable (RxSight)	RxLAL		Silicone 3-Piece (Bag)	NA / 1.43 +10.0 to +30.0	Adjustable diopteric power using UV light after implantation to correct for refractive misses Can also adjust astigmatism correction			
Akreos (B+L)	AO60 (Akreos AO) MI60P (Akreos MICS)	a de la compañía de	Hydrophilic Acrylic (Bag)	NEUTRAL/1.46 0.0 to +30.0 118.5/118.0   119.1/118.4 (AO) (MICS)	AO allows for 4-point scleral fixation; MICS can be implanted through small (1.8 mm) incisions for microincision cataract surgery			
AcrySof Multipiece (Alcon)	MA60AC (ant convex) MA50BM (post convex)		Acrylic 3-Piece (Sulcus)  ®	-0.20 / 1.55 +6.0 to +30.0 119.2/118.8   118.9	Small haptics of the 3-Piece IOLs minimize contact with uvea in the sulcus; MA60AC is anterior asymmetric biconvex, MA60BM posterior convex			
CZ70BD (Alcon)	CZ70BD		PMMA	-0.20 / 1.49 +10.0 to +30.0 118.8	Nonfoldable and rigid; good for open-sky cases with PKPs; eyelets for fixation			
SofPort (B+L)	LI61AO (preload) LI61SE (non-preload)		Silicone 3-Piece  ®	NEUTRAL/1.43 0.0 to +34.0 118.7/118.0	Can be placed in the sulcus or bag			
CT LUCIA (Zeiss)	CT Lucia 202 <sub>(spheric)</sub> CT Lucia 602 <sub>(aspheric)</sub>		Acrylic 3-Piece (Bag)	EITHER/ 1.49 +4.0 to +34.0 117.7	Very sturdy Polyvinylidene Fluoride (PVDF) 30G haptics; good for unstable bags and fixation			
			Extended De	pth of Focus				
Vivity (Alcon)	DFT015/DFTx CCWET <u>0/x</u> (Clareon)		Acrylic 1-Piece (Bag)  ©	-0.20 / 1.55 +10.0 to +30.0 119.2/118.8	Distance and intermediate vision; minimal near; non-diffractive "X-Wave technology"			
Symfony (Tecnis, J&J)	ZXR00/ZXR00V ZXT <i>x</i> /ZXW <i>x</i> (toric)		Acrylic 1-Piece (Bag)  ©	-0.27 / 1.47 +5.0 to +34.0 119.3/118.8	Distance/intermediate; blocks 360-460 nm; diffractive; "V" and "W" models are "OptiBlue"			
Synergy (Tecnis, J&J)	DFR00V DFW <i>x</i> (toric)	P	Acrylic 1-Piece (Bag)  ©	-0.27 / 1.47 +5.0 to +34.0 119.3/118.8	Has elements of EDOF (near and intermediate) and diffractive multifocal; with OptiBlue UV filter			
Eyhance (Tecnis, J&J)	DIBOO DIUx <sub>(toric)</sub>		Acrylic 1-Piece (Bag)  ©	-0.27 / 1.47 +5.0 to +34.0 119.3/118.8	FDA classified and priced as a <i>monofocal</i> with some intermediate using EDOF technology			
	1 .		Multi					
PanOptix (Alcon)	TFAT00/TFNT00 TFATx/TFNTx (toric) CCWTT <u>0/x</u> (Clarion)		Acrylic 1-Piece (Bag)  ©	-0.20 / 1.55 +6.0 to +30.0 119.1/118.7	Trifocal with true near; rings produce glare and halos around lights, especially at night			
Tecnis Multifocal	ZxB00 ZLUx (toric)		Acrylic 1-Piece (Bag)  ©	-0.27 / 1.47 +5.0 to +3.0 119.3/118.8	Bifocal with up to +4.00 D add; ZKB00, ZLB00, and ZMB00 for +2.75, +3.75, and +4.00 ADD			
Accommodating								
Crystalens; Trulign Toric (B+L)	AO1UV/AO2UV BL1UT <sub>(toric)</sub>		Silicone (Bag)	NEUTRAL/1.43 +4.0 to +33.0 119.1	Theoretically accommodates (vitreous pushes lens forward) to provide near vision; evidence for their efficacy was limited, they are not frequently used			

## (DTL) 119.1 efficiency was infinited, they are not freq \*Asphericity, Power Range, and Optical/US A-Constant; this list is not comprehensive and is not medical advice; it is for medical education reference; cartridges with the most compatibility by the Alcon Monarch system are encircled under "Ctg" Key: <u>Ux</u> = spherical/toric; example: CCWTT0<u>x</u> = CCWTT0 (monofocal) and CCWTTx (toric); © = preloaded or proprietary

Update Version: October 2022



## **Common Posterior Chamber IOLs**

## **Other FDA Approved and Currently Marketed\* IOLs**

	r			-	
Platform	Models	Picture	Material	Asph/n/D/A	Notes
Apthera IC-8 (AcuFocus)	Apthera IC-8	0	Acrylic 1-Piece (Bag)	Asph/1.483 +10.0 to +30.0	Optic has a nondiffractive 3.23-mm diameter opaque mask with a 1.36-mm central aperture to
				120.5/120.15	increase depth of field; aim -0.75
PreVue	PREVUEC (clear)		Acrylic 1-Piece	NEUTRAL/1.52	Haptics are modified-C loop, blue PMMA tips to
(B+L)	PREVUEY (yellow)		(Bag)	+6.0 to +30.0	increase surgeon visibility within the inserter and
(דר)		レレ		118.5/118.4	the eye; preloaded; minimizes "hand holding"
iPure (BVI)	B1PC (clear) B1PY (yellow) B3PC (clear, 3-piece)		Acrylic 1/3- Piece	-0.18	Very similar to PreVue (B+L) and like PreVue also
				+6.0 to +30.0	licensed from Hoya in the US. Unlike PreVue it has
			(Bag)	118.4	-0.18 of SA correction
iSert	iSERT 230 <sub>(clear)</sub> iSERT 231 <sub>(yellow)</sub>		Acrylic 3-Piece (Bag)	-0.18/1.52	Blue PMMA modified-C loop haptics, corrects the
(B+L)				+6.0 to +30.0 118.6/118.4	spherical aberration of the cornea, preloaded
Tecnis 3 Piece	ZA9003 (acrylic)		Acrylic/	-0.27 / 1.47	5°/10° of angulation, PMMA haptics
(Tecnis, J&J)	Z9002 (silicone)	$(\bigcirc)$	Silicone	+10/6 to +30.0	- ,,
(******	AAB00				3-piece option has 5° angulation and PMMA
Sensar	AR40e (3-Piece) AR40E (3-piece, low D) AR40M (3-piece, minus)		Acrylic (Bag)  © (Sulcus)  ®	Spherical/1.47 -10.0 to +30.0 119.0  118.7 (1P) (3P)	haptics, one of the rare PCIOLs with minus
(Tecnis, J&J)					powers, the AR40M ranges from -10.0D to +1.5D
					Applanation A-Constant is 118.4 for all Sensars
UltraSert	AU00T0 (yellow)		Acrylic 1-Piece	Negative/1.55	Preloaded version of AcrySof IQ
(Alcon)	ACUOTO (clear)		(Bag) ©	+6.0 to +30.0	
			(508)] ©	119.9 /118.7	
AcrySof		$\bigcap$	Acrylic 3-Piece	Negative/1.55	Meniscus, also a minus powered option
EXpand (Alcon)	MA60MA		(Sulcus) ®	-5.0 to +5.0 119.1/118.7	
	RAO200E (EMV)			Pos/Neut/Neg	Hydrophilic acrylic IOLs; EMV = enhanced
RayOne	RAO600C (aspheric) RAO100C (spheric)	SOJ	Acrylic 1-Piece (Bag)	+8.0 to +34.0 118.6	monovision like Eyhance but +SA
(Rayner)					n = 1.46 (ABBE 56); small injector nozzle
ReSTOR	(0)			Negative/1.55	Bifocal with +3.00 D add power, anterior
Multipiece	MN6AD1		Acrylic 3-Piece	+6.0 to +34.0	diffractive rings (most are posterior)
(Alcon)			(Bag)  ®	119.2	
	S122UV (12.50mm, 4.4°) L122UV (13.75mm, 3.7°)	SO.	PMMA 1-Piece (AC)	NEUTRAL +5.0 to +30.0 115.8/115.8	Lathe cut and polished from PPMA with an UV
B+L AC					absorbing compound added, S122UV for
(B+L)					10.00mm – 11.49mm WTW, L122UV for 11.50mm
					– 12.25mm WTW
Alcon AC (Alcon)	MTAxU0	NON	Silicone 3-Piece	-0.20 +5.0 to +30.0 <sup>115.3</sup>	MTA3U0 (12.5mm), MTA4U0 (13.0mm), and
					MTA5U0 (13.5mm) for different WTW
				Neutral	Hydrophilic acrylic, with a refractive segmental
SBL3	SBL3	87 077	Acrylic 1-Piece	+10.0 to +36.0	multifocal +3.00 ADD power
(LensTec)	-		(Bag)	118.0	
	MICL <i>x</i> TMICL <i>x</i> (toric) VICMO <i>x</i> (EVO)	$\sim$	Collamer (Sulcus)	Negative -3.0 to -18.0	ICL (implantable collamer lens) placed in sulcus
Vicion		• •			behind iris vaulted $1.0 \pm 0.5$ K thicknesses; the old
Visian, EVO Visian (STAAR)					version required LPIs, the new EVO version has a
	VICM5 <i>x</i> (EVO+)				360 $\mu$ m central port to eliminate this need
	VTICMOx (EVO toric) VTICM5x (EVO+ toric)				EVO+ has a 0.1 to 0.5 mm larger optic
					x represents vault sizing
Verisyse			PMMA	Negative	Attaches to the iris by "enclavation"; more visible
(I&I)	VRSM <i>x</i> US		(AC)	-5.0 to -20.0	than the ICL; $x = 5$ or 6 for mm sizing
. ,		<u> </u>	• •	116.8	

\*Based on current year product catalogs and spec sheets

- B 3-Piece IOLs, and high powers (35 D 40 D)
- © Any a©rylic
- O Most acrylics, some IOLs have a 25 D limit

